




WHITE PAPER

Optimizing Development and Scale-Up of Insoluble, Microbially Expressed Biologics

STEVE LOFTUS, PhD,
Microbial Business Steering Group,
FUJIFILM Biotechnologies

FUJIFILM
Value from Innovation



Inclusion bodies, historically considered undesirable byproducts of microbial expression systems, are now gaining attention for their potential to produce a range of difficult-to-express proteins in *E. coli* systems. The current biologics landscape is largely dominated by antibodies produced using CHO systems. Despite this, microbial systems, particularly *E. coli*, are still integral in enabling the production of more complex molecules less readily expressed in mammalian cells. *E. coli* can produce molecules that are soluble, expressed in the cytoplasm, fully folded, and active upon production.

However, *E. coli* can also produce inactive misfolded proteins that form masses of amorphous unfolded protein within the cell. While these misfolded proteins have long been considered an undesirable byproduct, more recent research has shown that these inclusion bodies may be valuable in increasing titers and simplifying primary separations.

Steve Loftus, PhD, Microbial Business Steering Group Lead for FUJIFILM Biotechnologies, explores the major challenges associated with the development of [microbial fermentation processes](#), as well as how our capabilities and expertise can help customers leverage inclusion bodies to improve their applications' CQAs.

Inclusion Bodies: Friend or Foe?

Over the past decade, data generated internally by FUJIFILM Biotechnologies is reflective of broader industry trends — namely, that the most commonly pursued route for *E. coli*-expressed biologics is soluble cytoplasmic products. The second most common selection is soluble periplasmic products, which accounted for just over a quarter of the products developed by FUJIFILM Biotechnologies; finally, approximately 18 percent of its *E. coli* products were those expressed insolubly in the cytoplasm. This disparity is chiefly due to the challenges of dealing with the amorphous masses produced by insoluble expression.

For non-mAb molecules on the market or in the late-stage pipeline today, many tend to more traditional biologics,

such as cytokines, hormones, and insulins. Yet newer biologic entities in earlier phases of development tend to be more complex molecules like antibody fragments, enzymes, or fusion proteins. This extra complexity can make it infeasible for these types of molecules to be expressed as a soluble protein in *E. coli*. As such, developing methods that enable the folding of these complex molecules is key to the success of these types of products. This is where a deeper evaluation of inclusion bodies comes into play.

The Advantages of Inclusion Bodies

Because they are produced as amorphous masses within the cell, inclusion bodies are highly resistant to proteolytic degradation, as the proteases cannot access the proteolysis sites within the molecule easily. These masses likewise protect the host cell from any proteins that may be inhibitory or toxic to the workings of that cell because they are not expressed in an active form. Moreover, applications can achieve extremely high initial titers for insoluble proteins because the folding machinery doesn't have to work — instead, ribosomes express a peptide chain repeatedly, building large amounts of protein within the cell due to the proteolytic resistance.

Additionally, inclusion bodies are generally stable across primary separations processes, whereas soluble proteins are more susceptible to being exposed to the contents of the cells or proteases, which can degrade the product if expressed in a soluble form during the primary



separation's stages. Employing insoluble expression can aid removal of many of the residual contaminants associated with microbial processes, such as endotoxins, lipids, or host cell proteins. Another key benefit is that upstream and downstream processes can be decoupled, allowing for the generation of substantial amounts of inclusion bodies, essentially in isolation from the rest of the process, which can be processed later when more downstream capacity is available.

The Challenges of Inclusion Bodies

The solubilization and refolding of inclusion bodies can be a complex process to develop, and the associated volumes can quickly become exceptionally large. To this end, it is important for process development scientists to consider the misfolding and aggregation of a product following the refolding process, as well as how to separate those misfolded and aggregated forms of the protein. These operations typically require more complex purification development and analytics to achieve their desired final purity.

No two refolds are likely to be the same, due to the individual folding kinetics of the proteins that are being expressed. While high initial titers have been an advantage of inclusion bodies, this also creates the potential for increased losses across the refold itself due to precipitation and misfolding. This can also cause losses during purification as these undesirable species need to be separated from the final product.

Pursuing Process Development for Insoluble Proteins

When it comes to both soluble and insoluble processes, FUJIFILM Biotechnologies has established a suite of technologies aimed at improving titers and addressing the key challenges that accompany these processes. For soluble expression, the FUJIFILM Biotechnology [Paveway™ PLUS](#) platform offers a flexible suite of modular workflows based on Paveway™, our proven technology platform for efficient microbial expression of proteins, using novel recombinant *E. coli* strains. This expression system offers industry-leading titers and a track record of more than 150 successfully expressed proteins. Development of a lead strain for recombinant protein expression can be achieved in as little as six weeks if selected on titer alone or 10 weeks if using high-throughput automated procedures to assess product quality. Paveway™ PLUS is particularly beneficial in scenarios where lead strain selection needs to ensure that no undesirable post-translational modifications have been introduced into a product because of the strain engineering performed.

After performing initial cloning and scale-down fermentation runs to produce up to eight top-performing strains, FUJIFILM Biotechnologies uses microscale harvest and primary separations to generate material for a high-level resin screen. This screen looks at basic conditions to determine if the product of interest will bind to one of two resins, enabling scale-up runs to generate enough material to

take forward into product analysis of these eight high titers producing strains. We evaluate intact mass, aggregation state, charge variance, and other key data points, and at the end of an 10-week workflow, present clients with the necessary data on each strain to inform decisions.

While this approach is highly effective for soluble molecules, FUJIFILM Biotechnologies has also developed a new high-throughput solution to slot into the Paveway™ PLUS workflow to enable work on products expressed insolubly. We have developed a two-stage workflow (Figure 1) based on liquid-handling robots, fed from lesser amounts of material

generated by Ambr® 250 fermenters, which provides enough data to perform lab scale refold experiments that provide material for later stages of the workflow.

The first stage is a solubilization study, which takes approximately five days and is split into a high-throughput automated stage and then the scale-up to produce material for stage two, which is the refold screen, followed by two rounds of Design of Experiments (DOE) screening and a subsequent lab scale-up approach to generate material for a resin screen.

Figure 1: Two-Stage Robotics Workflow




Table 1: Example Refold Buffer Screening Components

Buffer/pH	Stabilizers	Redox	Monovalent salt	Divalent salt/chelator
Buffer 1, pH 8.0	Stabilizer 1	Reductant 1	Salt 1	Divalent salt 1
Buffer 2, pH 6.5	Stabilizer 2	Redox Pair 2	Salt 2	Divalent salt 2
Buffer 3, pH 10.0	Stabilizer 3	Redox Pair 3	No monovalent salt	Divalent salt 3
	Stabilizer 4	No redox		Divalent salt 4
	Stabilizer 5			Chelator 1
	Stabilizer 6			No metal salt/chelator
	Stabilizer 7			
	Stabilizer 8			
	Stabilizer 9			
	No stabilizer			

The amount of data generated through these stages, prior to resin screening, affords FUJIFILM Biotechnologies considerable insight into whether the refold selected still produces the desired results at a larger scale, as well as the best buffer to utilize for subsequent stages. The number of factors evaluated in this process is extensive, as seen in Table 1, which outlines our refold screening evaluation of a “typical” refold buffer composition divided into five categorical factors — buffer/pH; stabilizers; redox; monovalent salt; and divalent salt/chelator.

In this example, FUJIFILM Biotechnologies analyzed three base buffers, nine different stabilizers, three redox pairs, two types of monovalent salt, four types of divalent salt, and a metal chelator. Overall, this equates to 2,160 possible combinations, compressed into 190 runs. Once complete, we can evaluate the best conditions identified during the runs, as shown in Figure 2.

We can also perform main effects modeling as part of this process, using results to build a predictive statistical model. Ultimately, final refold buffer selection is made by our DoE software, which uses the desirability function of its prediction profiler to select factors that maximize the refold index. This profiler gives components for a “best bet” buffer, which can be further investigated and optimized at lab scale. Overall, the enhanced insoluble product


workflow developed by FUJIFILM Biotechnologies adds approximately three weeks to Paveway™ PLUS’s existing 10-week timeline, enabling customers to access valuable data to inform later process development quickly.

Achieving Scale-Up for Insoluble Protein Expression

A typical downstream insoluble process starts with a large refold followed by many of the operations found in a soluble process, such as depth filtration to remove insoluble material, followed by a number of orthogonal chromatography steps often requiring the use of large-scale columns. While we currently possess a refold capacity of up to 10,000 liters utilizing traditional downstream approaches, it is moving toward more connected, continuous approaches to insoluble processing. Key to this transition is our proprietary, in-house developed SymphonX™ drug substance purification technology platform, a single multifunctional, multi-use bioprocessing system utilizing one disposable flow path capable of running filtration, chromatography, tangential flow filtration, and viral filtration. A one-stop shop for downstream processing, the FUJIFILM Biotechnologies SymphonX™ system can also perform advanced buffer management, including inline dilution and conditioning;

Figure 2: Scatter Plot Matrix of Stage 1 Refold Results

This scatter plot shows the results of a series of Design of Experiments (DOE) runs, organized by buffer component, indicating a clear positive effect for refold buffers.



moreover, because it is fully automated, it can be connected to other SymphonX™ systems to establish a full-scale continuous process.

FUJIFILM Biotechnologies has worked to streamline its process development and optimize scale-up for insoluble expression, focusing on tweaking parameters for different expression routes, optimizing continuous centrifugation at the pilot scale, and harmonizing process development equipment for use in later manufacturing. It has performed similar optimization across the process's harvest, wash, solubilization, and refolding steps, screening for optimal conditions for washing, for example, or employing advanced analytics to assess refolding properties. By prioritizing high-throughput technology, connected, continuous processing, and seamless scaling, we are working to establish a workflow for insoluble proteins that offers a complete solution to surmounting the challenges associated with insolubly expressed biologics.

About the Author



STEVE LOFTUS, PhD

Microbial Business Steering Group,
FUJIFILM Biotechnologies

Dr. Loftus leads the Business Steering Group for Microbial Services at FUJIFILM Biotechnologies, a group

that looks to define the strategic growth strategy of the offering. He has over 15 years of experience in the development and manufacture of microbial-based biologics at FUJIFILM Biotechnologies. He has a doctorate in Biochemistry and Biophysics from the University of York.

Partners for *Life*

"Partners for Life" represents a transformative approach to development and manufacturing, emphasizing relationships founded on trust and transparency — founded in people-centric values, transformative science and innovation, and unprecedented delivery.

CDMO
LEADERSHIP
AWARDS 2024

**Contact us to discuss
your science.**

USA

College Station, Texas
3940 Fujifilm Way
College Station, TX 77845
+1 979 431 3500

USA

Holly Springs, North Carolina
100 Biotechnology Avenue
Holly Springs, NC 27540
(Coming online 2025)

USA

**Research Triangle Park,
North Carolina**
101 J Morris Commons Lane
Morrisville, NC 27560
+1 919 337 4400

USA

Thousand Oaks, California
2430 Conejo Spectrum Street
Thousand Oaks, CA 91320
+1 805 699 5579

DENMARK

Hillerød
Biotek Allé 1
3400 Hillerød
+45 7741 6000

UNITED KINGDOM

Teesside
Belasis Avenue
Billingham, TS23 1LH
+44 1642 363511

JAPAN

Tokyo
R7 Building, 7-12-2 Roppongi
Minato-ku, Tokyo 106-0032
+81 3 6871 7736



FUJIFILM Biotechnologies

fujifilmbiotechnologies.fujifilm.com

Email: contactfdb@fujifilm.com



MKT-0179-WP-Rev01
© 2025 FUJIFILM Diosynth Biotechnologies.
All Rights Reserved.